## **UNIVERSITY OF MISSOURI** Y-90 SIR-TeX MICROSPHERE ADMINISTRATION FORM

Patient sticker here To include:		PATIENT IDENTIFICATION VERIFIED by (2 required):  □ Name □ Hospital I.D. □ Date of Birth □ Personal Recognition □ Photo I.D.
Patient is:  Male Female, declared not currently pregnant or breastfeeding  RADIOPHARMACEUTI		VERIFY PATIENT IS NOT PREGNANT by:  □ Declares menstruating, not pregnant, or to be postmenopausal <u>OR</u> □ Negative pregnancy Test <u>OR</u> □ Rendered sterile, i.e. Hysterectomy, Tubal Ligation, etc., <u>OR</u> □ Verify patients' age is ≤ 8 or ≥ 60 years    CAL PRESCRIPTION:
(To be completed by the Authorized User)		
PRE-IMPLANT DOCUMENTATION:		RADIOPHARMACEUTICAL DOSAGE VERIFICATION:
Treatment site: □ Right lobe □ Left Lobe □ Entire liver		Lot No
Date and Time of Procedure:		Prepared Activity [KK]:mCi
Radionuclide/Physical Form: <u>Y-90 SIR-Spheres resin microspheres</u>		Date and time:
Manufacturer: Sirtex Wilmington, LLC		Is the MUHC SIR-Spheres Y-90 Resin Microspheres Verification
Route of Administration: <u>Intra-arterial</u>		Worksheet attached? □ Yes □ No
Prescribed Activity:mCi		
OR activity delivered at stasis		Pre-/Post-Administration checks performed by signature:
% Lung Shunting: %		Date:
Maximum Shunting Activity:mCi  Is the SIR-Spheres Activity Calculator (SMAC) attached? □ Yes □ No Authorized User Signature:		Note: SirTex Sir-Spheres package insert contraindications advise a maximum permissible dose limit to the lung due to arterial shunting of 30 Gy.
Date:		
POST-IMPLANT:		
		and dated by an AU within 24 hours.
		Date and Time of Delivery:
Calculated Activity Outside of the Liver from Shunting: mCi Dose Delivered Due to Shunting: Gy (<30Gy)		
Stasis? □ Yes □ No Emerging patient conditions? □ Yes □ No  Provide explanation:		
rrovide explanation:		
Authorized User Signature:		
Date:		
PATIENT RELEASE:		
□ Patient reading at 1 m is <2.0mR/hr □ Post-implant instructions given		
Final review by NM (Signature):		Date:/

RETAIN THIS RECORD FOR 3 YEARS: REQIORED BY 10 CFR 35.2040